

**Baxter**



Introducing the New

# Welch Allyn Connex 360 Vital Signs Monitor



It's About **Time**

# It's About Time for Future-Ready Vitals Management

## Meet the **Welch Allyn Connex 360** Vital Signs Monitor

With an upgradeable and expandable design, the **Connex 360** features optimized clinician workflows for connected monitoring of adult, pediatric and neonate patients\*.

- 12" color touchscreen
- Continuously measure SpO<sub>2</sub>, respiration rate\*, and pulse rate with refreshed data once per second
- Customizable intervals monitoring programs for timed NIBP measurements and BP/PR averaging
- Take and document a full set of vitals in less than one minute
- Profile-free intuitive navigation
- Apply customizable Early Warning Score protocols (eg. MEWS, PEWS, NEWS2) to help detect patient deterioration and provide clinical direction
- Select Masimo or Nellcor SpO<sub>2</sub> technologies
- Document up to 20 fields of patient information for more comprehensive bedside workflows

\*Digitally measure respiration rate with **Masimo RRRp**, for patients greater than 2 years old.

## Designed for Optimized Workflows

- 1 Our largest and most responsive color touchscreen
- 2 Imprivata badge access for streamlined, secure clinician authentication with single sign-on
- 3 Flexible monitoring capability with easy-to-use, intuitive navigation between spot check, timed intervals monitoring and averaging programs



## Help Improve Patient Safety



Let's help give clinicians time back in their day and help improve patient outcomes with EMR connectivity.

Connected vitals solutions give clinicians access to patient information to help inform clinical decisions and reduce transcription errors.



## DESIGNED WITH A COMPREHENSIVE SECURITY FOCUS



The new **Connex 360** platform and **DeviceBridge** solution were developed with security in mind from conception to deployment. Security-forward features help ensure protection and encryption of data at rest on device and in transit. Enhanced security features include:

- Trusted Platform Module (TPM) to help enhance security of the device hardware and software
- Certificate-based mutual authentication using Baxter PKI or customer provided certificates
- FIPS-approved encryption algorithms of transmitted data
- Supports TLS 1.3 and WPA3 Security Protocols

## DEVICES

Part Number	Description
94CXT	Welch Allyn Connex 360, Nellcor SpO <sub>2</sub> , SureTemp
94MXT	Welch Allyn Connex 360, Masimo SpO <sub>2</sub> , SureTemp
94RXT	Welch Allyn Connex 360, Masimo RRp SpO <sub>2</sub> , SureTemp
95CXT	Welch Allyn Connex 360, Wi-Fi, Nellcor SpO <sub>2</sub> , SureTemp
95MXT	Welch Allyn Connex 360, Wi-Fi, Masimo SpO <sub>2</sub> , SureTemp
95RXT	Welch Allyn Connex 360, Wi-Fi, Masimo RRp SpO <sub>2</sub> , SureTemp

## KITS

C360KIT95CXTBMBJX	C360 Kit: 95CXT-B, 9000-C360 Mobile Stand, JADAK 7000-916HS 1RS Scanner
C360KIT95MXTBMBJX	C360 Kit: 95MXT-B, 9000-C360 Mobile Stand, JADAK 7000-916HS 1RS Scanner
C360KIT95RXTBMBJX	C360 Kit: 95RXT-B, 9000-C360-Mobile Stand, JADAK 7000-916HS 1RS Scanner
C360KIT95CXTBWLJBP	C360 Kit: 95CXT-B, 9000-GCX Wall Mount, JADAK 7000-916HS 1RS Scanner, FlexiPort 10-foot Fast BP Hose 4500-35
C360KIT95MXTBWLJBP	C360 Kit: 95MXT-B, 9000-GCX Wall Mount, JADAK 7000-916HS 1RS Scanner, FlexiPort 10-foot Fast BP Hose 4500-35
C360KIT95RXTBWLJBP	C360 Kit: 95RXT-B, 9000-GCX Wall Mount, JADAK 7000-916HS 1RS Scanner, FlexiPort 10-foot Fast BP Hose 4500-35

## MOUNTING OPTIONS

9000-GCX	GCX VESA Wall Mount
9000-C360	Welch Allyn Connex 360 Mobile Stand

## SCANNERS

7000-916HS	JADAK Barcode Scanner and Holster Kit
7000-916HSR	JADAK Barcode/High Frequency RFID Scanner and Holster Kit
7000-916HS1RS	JADAK Barcode/Dual Frequency RFID Scanner and Holster Kit
7000-CS6080	Zebra Barcode Scanner and Holster Kit

Full line of Accessories available. See website for details.



With Imprivata Single Sign-On, clinicians can log into the device with a simple tap, no passwords needed.



Let's get your vitals management future-ready.

Contact your Baxter account representative today for a demo or visit [hillrom.com/connex360](https://hillrom.com/connex360).

**Rx Only.** For safe and proper use of the product mentioned herein, please refer to the appropriate Operator's Manual or Instructions for Use.

### Baxter.com

Baxter, Connex 360, SureTemp and Welch Allyn are trademarks of Baxter International Inc. or its subsidiaries.

Masimo and RRp are trademarks of Masimo Corporation.

Nellcor is a trademark of Covidien LP.

Wi-Fi is a trademark of Wi-Fi Alliance.

Any other trademarks, product names or brand images appearing herein are the property of their respective owners.

US-FLC215-250001 (v3.0) 12/2025