


BLOOD DRAWING CHAIRS



PADED BLOOD DRAWING CHAIR

- Overall Dimensions: 36" H X 30.5" W X 28" D
- Overall Weight: 35 lb
- Weight Capacity: 400 lb

- Padded for optimal patient comfort
- Armrest that can lock into the correct position
- Easy-to-clean antimicrobial polyurethane leather
- Welded powder-coated frame with reinforced stress points



**ARMREST
HEIGHT ADJUSTMENT:**
24.4" – 34.2"

BLACK
ADI997-01-BLK

MIDNIGHT EXP.
ADI997-01-BLU

500 LB CAPACITY

Optional Accessory

DRAWER

for Padded
Blood Drawing Chair

BLACK # ADI997-01-D-BLK **MID. EXP.** # ADI997-01-D-BLU



PADED BLOOD DRAWING CHAIR

With Drawer

- Amply sized 4" x 6" drawer for various medical supplies

BLACK
ADI997-01D-BLK

MIDNIGHT EXPRESS
ADI997-01D-BLU

BARIATRIC PADDED BLOOD DRAWING CHAIR

- Overall Dimensions: 40" H X 35" W X 30" D
- Overall Weight: 48.5 lb

- Securely holds up to 500 lb
- Padded for optimal patient comfort
- Armrests that can lock into the correct position
- Easy-to-clean antimicrobial polyurethane leather
- Welded powder-coated frame with reinforced stress points



**ARMREST
HEIGHT ADJUSTMENT:**
27.8" - 38.5"



BLACK

ADI997-01-BLK-HD



MIDNIGHT EXP.

ADI997-01-BLU-HD

500 LB CAPACITY

Optional Accessory

DRAWER

for Bariatric Padded
Blood Drawing Chair

BLACK

ADI997-01-D-BLK

MID. EXP.

ADI997-01-D-BLU



BARIATRIC PADDED BLOOD DRAWING CHAIR

With Drawer

- Amply sized 4" x 6" drawer for various medical supplies



BLACK

ADI997-01D-BLK-HD




MIDNIGHT EXPRESS

ADI997-01D-BLU-HD

LUXE UPHOLSTERED BLOOD DRAWING CHAIR

- Padded for optimal patient comfort
- Armrest that can lock into the correct position
- Easy-to-clean antimicrobial polyurethane leather
- Welded powder-coated frame with reinforced stress points
- Overall Dimensions: 35" H X 28" W X 25" D
- Overall Weight: 54 lb



BLACK

ADI997-02-BLK



MIDNIGHT EXPRESS

ADI997-02-BLU



FOLDING PADDED BLOOD DRAWING CHAIR

ADI997-03-BLK

- Can be effortlessly transported
- Foldable design to facilitate storage
- Easy-to-clean antimicrobial upholstery
- Seat designed with a thick cushion for patient comfort
- Equipped with nylon glides on the legs to avoid scratching your floors
- Overall Dimensions: 37" H X 25.5" W X 24" D
- Overall Weight: 28 lb



EASILY FOLDABLE

