

NOBLE-LINE TREATMENT TABLES



NOBLE-LINE ADJUSTABLE TREATMENT TABLE

With Paper Dispenser

- Foldable leg design helps provides space-saving storage
- Stainless-steel frame for longevity and easy cleaning
- Includes paper roll dispenser and paper cutter
- Overall Dimensions: 31.5" H X 27.5" W X 75.5" L
- Total Weight: 60 lb
- Adjustable backrest
- Rubber feet help prevent floor scratches
- Antimicrobial mattress provides a hygienic environment



BLACK

ADI996-01-BLK

MIDNIGHT EXPRESS

ADI996-01-BLU



SANITATION

Antimicrobial mattress provides safe and hygienic environment, and reduces the risk of cross-contamination

2-INCH PADDED
FOAM PU LEATHER
UPHOLSTERY



FOLDABLE LEGS


KIT #	DESCRIPTION	INCLUDED SKU _s	COLOR	
ADI996-01-DR3-BLK	Adjustable Exam Table + Three Drawer Set	ADI996-01-BLK ADI996-01-DR-3	Black	
ADI996-01-DR3-BLU	Adjustable Exam Table + Three Drawer Set	ADI996-01-BLU ADI996-01-DR-3	Blue	

Includes
PAPER CUTTER

Includes
PAPER DISPENSER

**ADJUSTABLE
BACKREST**

Can be positioned at
different heights for the
patient's optimal
comfort



Also available with
3 DRAWER SET

ADI996-01-DR-3

NOBLE-LINE ADJUSTABLE EXAM TABLE


With Shelf and Exam Paper Roll Dispenser

- Foldable leg design helps provides space-saving storage
- Stainless-steel frame for longevity and easy cleaning
- Includes paper roll dispenser
- Overall Dimensions: 31.5" H X 27.5" W X 74.8" L
- Total Weight: 73.6 lb
- Adjustable backrest
- Rubber feet help prevent floor scratches
- Antimicrobial mattress provides a hygienic environment
- Backrest adjust to angles: 0°, 32.5°, 45°, 54.5°, 62°, and 67°

2-INCH PADDED
FOAM PU LEATHER
UPHOLSTERY

ADI996-02-BLK

LB
440 LB
CAPACITY



PAPER DISPENSER



ADJUSTABLE BACKREST



FOLDABLE LEGS



SHELF

NOBLE-LINE MAHOGANY WOODEN EXAM TABLE

With Full Shelf

- Equipped shelf provides storage space for various medical equipment
- Wooden frame for longevity and aesthetic appeal
- Integrated paper roll dispenser
- Overall Dimensions: 31.5" H X 28.5" W X 75.5" L


ADI996-03-MA

- Total Weight: 112 lb
- Adjustable backrest
- 2-inch padded mattress delivers optimal comfort



CONVENIENT

Paper roll dispenser will help reduce the risk of cross-contamination between patients



**BUILT-IN SHELF PROVIDES
AMPLE STORAGE**

NOBLE-LINE BLACK STRAIGHT LINE TREATMENT TABLE

- Wooden frame for longevity and aesthetic appeal
- Rounded corner top ensures patients' safety
- 2.25-inch padded mattress delivers optimal comfort
- Overall Dimensions: 31" H X 24" W X 72" L
- Total Weight: 77 lb


ADI996-04-BLK

2.25-INCH PADDED
FOAM UPHOLSTERY



PROTECTION

Rounded corner top to ensure safety and long-term use



DURABLE

Reinforced joint corners enhance stability and longevity