

SHARPS & NEEDLE DISPOSAL CONTAINERS



SHARPS AND NEEDLE DISPOSAL CONTAINER

1 Quart

- Dimensions Per Container: 6" H X 4.3" W X 4" D
- Weight Per Container: 0.19 lb

- Designed to be puncture resistant
- Made from robust, heat-resistant PP material
- Features a split-lid design and dual openings
- Temporary closure and snap locks provide ultimate protection
- Crafted with a transparent lid to see the filling progress

1-Pack

ADI998-01-01



3-Pack

ADI998-01-03



6-Pack

ADI998-01-06



SPLIT-LID DESIGN



SHARPS CONTAINER WITH MAILBOX-STYLE HORIZONTAL LID


5 Quarts

- Dimensions Per Container: 11" H X 10.6" W X 4.5" D
- Weight Per Container: 0.66 lb

- Designed to be puncture resistant
- Made from robust, heat-resistant PP material
- Mailbox-style lid allows the disposal of needles without direct contact
- Crafted with a transparent lid to see the filling progress
- Outfitted with locks to ensure safe transportation and sharps disposal



LID OPENING



1-Pack

ADI998-02-01



2-Pack

ADI998-02-02



SHARPS CONTAINER WITH FLIP-OPEN LID

8 Quarts

- Dimensions Per Container: 9.8" H X 8.2" W X 7" D
- Weight Per Container: 0.73 lb

- Features a flip-open lid design
- Designed to be puncture resistant
- Made from robust, heat-resistant PP material
- Crafted with a transparent lid to see the filling progress
- Accommodates larger syringes and has a significant holding capacity

1-Pack

ADI998-03-01



2-Pack

ADI998-03-02



FLIP-OPEN LID DESIGN



SHARPS SHUTTLE CONTAINER

With Locking Mechanism

- Dimensions Per Container: 6.5" H X 1.2" D
- Total Weight: 0.26 lb

- Designed to be puncture resistant
- Made from robust, heat-resistant PP material
- Allows for one-handed operation
- Cap is specifically designed to be leak-resistant, providing secure containment
- Convenient size easily fits into bags for EMS, Fire and Rescue, and Nursing

6-Pack

ADI998-04-06



LEAK-RESISTANT CAP



SHARPS CONTAINER ROUND SHAPED

1.5 Quarts

- Dimensions Per Container: 6" H X 5.5" W X 4.3" D
- Weight Per Container: 0.25 lb

1-Pack

ADI998-05-01



2-Pack

ADI998-05-02



SAFETY



Outfitted cover ensures that the contents stay safely inside the container