

OPHTHALMIC CATALOG



WelchAllyn®

Advancing Frontline Care™

Welch Allyn Advancing Frontline Care

In ophthalmic practices, the people you see every day are more than just your patients—they are your friends and neighbors. And when patients visit your office, they're not just going to the doctor—they're going to their doctor.

For nearly 100 years, Welch Allyn has been right there with you, providing innovative, affordable healthcare products and solutions that allow you to more effectively diagnose and treat your patients—whether you are caring for infants, the elderly, or any age in between.

At Welch Allyn, "innovation" means providing you with state-of-the-art products designed to meet your need for simplicity and efficiency. So you can see more patients, perform more procedures, and provide more complete on-site care.



OPHTHALMIC

Introduction	3
PanOptic™ Ophthalmoscope	5
Aperture and Filter Guide	4
3.5 V Coaxial Ophthalmoscopes	6
LED Replacement Lamps	7
Retinoscopes	8
PocketScope™ and CompacSet™	9
Diagnostic Systems and Sets	10
767 Wall Transformer	13
Binocular Indirect Ophthalmoscope	14
SureSight™ Autorefractor	15
Patient Monitors and Systems	16
Vital Signs Assessment	16, 20
Surgery and Recovery Equipment	18
Medical Exam Lighting	22
Replacement Lamps	24
Replacement Batteries	26
Welch Allyn Locations	27

For more information on the products in this catalog, visit our website at www.welchallyn.com.

A Tradition of Visionary Thinking

Since our inception in 1915, Welch Allyn has focused on developing instruments for busy clinicians to diagnose and treat their patients. This has allowed us to build a strong reputation as a leading manufacturer of medical diagnostic instruments. In fact, Welch Allyn has created innovations that have set the standards by which all others are compared. Here are just a few important milestones we've achieved.

- 1915 World's first hand-held, directly illuminated ophthalmoscope
- 1973 First use of halogen illumination in a medical instrument
- 1983 First coaxial optical system in a direct ophthalmoscope
- 1989 First 68-lens automatic diopter progression and display ophthalmoscope
- 1999 First lithium ion battery technology in hand-held ophthalmic instruments
- 2001 The PanOptic™ Ophthalmoscope, a revolutionary instrument that provides a 5X larger view of the fundus in an undilated pupil
- 2004 The Elite™ Retinoscope provides 2.8X the brightness compared to our standard retinoscopes
- 2010 LED Replacement Lamps provide a whiter, brighter light than halogen with 10,000 hours of lamp life

The advances listed here not only represent significant advances in technology, but also are durable, easy-to-use products that will save you time on the front lines of care. Our focus has always been on making better products to save you time, so you can focus on what's really important—your patients.

In this catalog you will find information on our line of ophthalmic diagnostic instruments. It includes information on product features and benefits, as well as model numbers to make ordering easy. It even includes information on replacement lamps and batteries. All Welch Allyn ophthalmic instruments are guaranteed by Welch Allyn against manufacturing defects. Welch Allyn will repair or replace, free of charge, any parts of its own manufacture proven to be defective through causes other than misuse, neglect, damage in shipment, or normal wear.

If you would like to place an order or have any questions on a product listed here, call 1.800.535.6663 to speak with a Welch Allyn Customer Service representative.

Visit us at www.welchallyn.com



PanOptic™ Ophthalmoscope

What makes the PanOptic™ Ophthalmoscope so revolutionary? Quite simply, the Welch Allyn Axial PointSource Optics. By converging the illumination to a point at the cornea, Axial PointSource Optics allow easy entry into even the smallest pupils. They then enable the illumination pathway to diverge to the retina, creating a very wide illuminated area of the fundus.

Simultaneously, the scope's unique viewing system allows you to see the illuminated area on virtually the same axis, creating one of the widest fields of view attainable in undilated ophthalmoscopy. This means you are actually viewing the fundus through the mirror that directs illumination into the eye, maximizing the alignment of the viewing axis with the illuminated area. The bottom line? The Welch Allyn PanOptic Ophthalmoscope provides a dramatic difference you can see right away. And that translates into easier, more complete fundoscopic exams for you—the first time and every time.

#11820 PanOptic Ophthalmoscope with Cobalt-Blue Filter and Corneal Viewing Lens (13.3 diopters)

#11870 Patient Eyecups (package of 5 reusable cups)

#03800-U Halogen HPX™ Replacement Lamp

#03800-LED LED Replacement Lamp

#11875 Corneal Viewing Lens (for use with 11820)

#05828-U Diagnostic Set Carrying Case (hard)

#41820 Apertures and Filters: 

ADVANTAGES OF USING THE PANOPTIC

- A 5X larger view of the fundus vs. standard ophthalmoscopes in an undilated eye
- Fast, easy entry into small, undilated pupils
- A staggering 25° field of view vs. the 5° field of view of standard ophthalmoscopes
- Greater working distance between practitioner and patient, improving comfort for all
- A 26% increase in magnification, making it easier to see retinal details

PATIENT EYECUP:

The eyecup establishes and maintains proper viewing distance, provides image stability, serves as a pivoting point for panning around the retina, helps with orientation, and screens out ambient light. The eyecup can be cleaned and is reusable

AXIAL POINTSOURCE OPTICS:

Easy entry into undilated pupils, resulting in a view of the fundus that's 5X larger than standard ophthalmoscopes in an undilated pupil

PATENTED GLARE EXTINGUISHMENT SYSTEM:

Prevents interference from unwanted glare and reflections

DYNAMIC FOCUSING WHEEL:


Continuous, smooth action for more precise control and optimum view (+20 to -20 diopter range)

PanOptic Ophthalmoscope
with Lithium-Ion Power Handle

APERTURE AND FILTER DIAL:

Allows you to select micro, small, and large spot sizes. Also includes a slit aperture and both a cobalt blue and red-free filter






WHAT YOU WILL SEE

The entire optic disc and surrounding vessels should be visible.

You should be realizing 5X the field of view you would see with a standard ophthalmoscope in an undilated eye.

This wide view makes it much easier to see important fundus structures.



APERTURE AND FILTER GUIDE FOR ALL WELCH ALLYN OPHTHALMOSCOPES :



MICRO-SPOT APERTURE: allows quick entry into small, undilated pupils.



SMALL APERTURE: provides excellent view of fundus through an undilated pupil.



LARGE APERTURE: standard aperture for dilated pupil and general examination of the eye.



FIXATION APERTURE: has graduated cross-hairs for measuring eccentric fixation or locating lesions.



SLIT APERTURE: used to determine levels of lesions and tumors.



COBALT FILTER: used with fluorescein dye to view small lesions, abrasions, and foreign objects on the cornea.



RED-FREE FILTER: Excludes red rays from examination field for easy identification of veins, arteries, and nerve fibers.



FILTER SWITCH: Rotate between red-free, polarized and unpolarized views. Use with the first 6 apertures shown yields a total of 18 possible apertures.

Coaxial Ophthalmoscopes

The Welch Allyn patented Coaxial Vision System makes ophthalmoscopy easier than standard aspheric ophthalmoscopes. With its combination of patented Coaxial optics, a unique new selection of filters and precise manufacturing techniques, the Coaxial Ophthalmoscope optics produce a shadow-free spot, easier entry into undilated pupils and a larger field of view.

Apertures and Filters: 

#11735
PATIENT SIDE



Prestige™ Coaxial-Plus Ophthalmoscope

- 2 filter selections provide optimum light levels for entering healthy eyes and those with cataracts
- 68 focusing lenses range from -30 to +38 diopters for better resolution
- Neutral filter provides the ideal amount of color-accurate light for entering healthy eyes without unwanted reflex
- Opacity setting provides the brightest light for easy entry

#11735 Prestige™ Coaxial-Plus Ophthalmoscope

#04900-U Halogen HPX Replacement Lamp

#04900-LED LED Replacement Lamp

#11730
PATIENT SIDE



AutoStep® Coaxial Ophthalmoscope

- 18 aperture/filter combinations for greater versatility
- 68 focusing lenses in single-diopter steps for precise resolution range of -30 to +38 diopters
- Sealed optics keep out dust and dirt for years of maintenance-free operation
- Rubber brow rest prevents scratching of eyeglasses

#11730 3.5 V Coaxial Ophthalmoscope

#04900-U Halogen HPX Replacement Lamp

#04900-LED LED Replacement Lamp

#11720
PATIENT SIDE



Coaxial Ophthalmoscope

- 18 aperture/filter combinations for greater versatility
- 28 focusing lenses with a range of -25 to +40 diopters
- Detect corneal abrasions with cobalt blue filter
- Red-free filter may be used with any aperture
- Polarizing filter virtually eliminates corneal reflection
- Sealed optics keep out dust and dirt for years of maintenance-free operation

#11720 3.5 V Coaxial Ophthalmoscope

#04900-U Halogen HPX Replacement Lamp

#04900-LED LED Replacement Lamp

COAXIAL OPTICAL SYSTEM:

Precise lamp filament and optical alignment enable maximum coaxial viewing for enhanced visibility and uniform illumination

IMPROVED LENS SELECTION DISC:

Offers 68 lenses with smooth, simple movement, and instant readout of actual diopter level

APERTURE SELECTION DIAL:

Offers five standard apertures

OPACITY SETTING:

Provides the brightest light for easy entry into eyes with cataract or other media opacity

NEUTRAL FILTER:

Provides the ideal amount of light for entering healthy eyes without unwanted reflex

3.5 V Coaxial Ophthalmoscope with Lithium Ion Power Handle





LED Replacement Lamps

30X the lamp life of a traditional halogen lamp in just one LED lamp

The Welch Allyn LED lamp provides up to 80% cost savings over the life of the lamp versus our equivalent halogen lamp, plus:

- Triple the examination time on a fully charged battery handle
- Consistent color even when light intensity is reduced
- LED color and patented circuitry designed specifically for intended examinations
- Lower heat build up for increased clinician and patient comfort
- 7-year guarantee

LED Lamp Part #	Device Part #	Description
 04900-LED	Traditional Ophthalmoscopes	
	11720	3.5 V Coaxial Ophthalmoscope
	11730	AutoStep® Coaxial Ophthalmoscope
	11735	Prestige™ Coaxial-Plus Ophthalmoscope
 03800-LED	PanOptic™ Ophthalmoscopes	
	11810	PanOptic™ Ophthalmoscope
	11820	PanOptic™ Ophthalmoscope/Blue Filter/Lens

LED Color Temperature Spectrum

Welch Allyn Ophthalmoscope—LED
3500 °K



Welch Allyn Ophthalmoscope—Halogen
3200 °K – 2800 °K



Dimming Cycle

LED




Consistent color output through dimming cycle

@100%

@25%

Halogen



Elite™ Retinoscope

SUPERIOR ILLUMINATION FOR FASTER AND EASIER REFRACTIONS

The Elite Retinoscope has an enhanced light output that makes refractions easier than ever by producing one of the brightest, sharpest, and cleanest streaks available. Features such as a convenient external focusing sleeve, sealed optics, and an optional polarizing filter make this one of the most user-friendly and innovative retinoscopes offered.

SUPERIOR STREAK QUALITY:

Enables easier and quicker refractions—even in patients with small pupils

CROSSED-LINEAR POLARIZING FILTER:

Eliminates bothersome glare from trial lenses, helping you perform a more precise exam at the flick of a switch

EXTERNAL FOCUSING SLEEVE:

The continuously rotating external focusing sleeve improves comfort and maneuverability no matter what size your hand is or how you prefer to manipulate the sleeve


DUST-FREE DESIGN:

Completely sealed optics keep the lenses cleaner for more effective examinations over a longer period of time

BRIGHTNESS VS. OUR STANDARD RETINOSCOPE:

ELITE OPTICAL SYSTEM:

The Elite Optical System provides 2.8X the brightness to the user and superior streak quality versus our standard retinoscopes



- #18245 Welch Allyn Elite Streak Retinoscope
- #08200-U Halogen HPX Streak Replacement Lamp
- #08300-U Halogen HPX Spot Replacement Lamp
- #18250 Set of four fixation cards quickly and easily attach to magnetic logo on instrument head for dynamic retinoscopy.

SPOT RETINOSCOPE INTERCHANGEABILITY

Simply insert a spot filament bulb in place of the streak filament bulb and the Welch Allyn Elite Retinoscope is converted into a spot retinoscope. This allows you to perform spot retinoscopy without having to buy a second scope.



Elite™ Retinoscope
with 719 Lithium-Ion Handle

Ophthalmic PocketScope™

Welch Allyn 2.5 V Ophthalmic PocketScopes are lightweight, compact diagnostic instruments for those who want high performance in a slightly smaller package. Welch Allyn PocketScopes use halogen illumination and have many other features found on 3.5 V sets. They are smaller in size but not in utility.



#12820 PocketScope Ophthalmoscope

#16220 PocketScope Retinoscope

Ophthalmic CompacSet™

This portable ophthalmic instrument is specifically designed to make examining patients easier. The Welch Allyn CompacSet features a retinoscope and an ophthalmoscope built into one compact unit.



#93300 Ophthalmic CompacSet with AA batteries and case



POCKETSCOPE SET

#94830 (94805* export) Ophthalmic PocketScope Set with rechargeable 2.5 V batteries, charging stand, and case



INDIVIDUAL POCKETSCOPE

- #12820 PocketScope Ophthalmoscope with AA batteries
- #16220 PocketScope Retinoscope with AA batteries
- #14220 PocketScope Transilluminator with AA batteries

ACCESSORIES

- #03900-U PocketScope Ophthalmoscope lamp
- #03400-U PocketScope Transilluminator lamp
- #08400-U PocketScope Retinoscope lamp
- #72600 Rechargeable 2.5 V PocketScope NiCad battery
- #72610 Rechargeable 2.5 V Ophthalmic CompacSet NiCad battery

*Available for export. Please specify country and voltage. Export ordering example: #94805 – Germany – 220 V

3.5 V Lithium Ion Rechargeable Handle

Twice the power, half the weight.

- Smaller diameter and half the weight of a full-size 3.5 V NiCad handle
- 120 minutes of operation between charges (compared to 60 minutes for standard NiCad handles)
- Charges in Universal Desk Charger or with removable AC adapter
- Compatible with all Welch Allyn 3.5 V instrument heads
- Low-battery indicator light

#71900 (#71905* export) Lithium-Ion Power Source with AC Adapter

#71930 (#71935* export) Lithium-Ion Power Source with Well Adapter

#71960 Lithium-Ion Battery



SOFT GRIP RHEOSTAT:
For light-intensity control

TEXTURED HANDLE:
Made of durable plastic-resin
blend construction



#71640 Desk Charger with
Lithium-Ion Handles
(heads sold separately)



#71980 Lithium-Ion
Well Adapter

Allows use of Lithium-Ion Power
Handle with traditional charging stand.



#71950 Removable
Pull-out AC
Charging
Module



Retractable AC charging module easily
plugs into any outlet.

RECYCLABLE BATTERY:
Please return expended
rechargeable batteries to
Welch Allyn for recycling.



*Available for export. Please specify country and voltage. Export ordering example: #71905 – Denmark – 220 V

3.5 V Nickel Cadmium Rechargeable Handles

A choice of convenient rechargeable power sources that fit all our 3.5 V instrument heads.

- Smooth, nonabrasive knurled finish ensures a no-slip grip
- Chrome-plated brass construction for durability and reliable long life
- 60 minutes of operation between charges
- Rechargeable and can also be converted for use with standard C-cell batteries
- Charges in Universal Desk Charger and numerous stand consoles and well chargers

#71670 3.5 V Well Rechargeable Handle

#71630 Universal Desk Charger and Two NiCad Handles (heads sold separately)

#71000-C 3.5 V Rechargeable/Convertible Handle (wall plug-in)

#71000-A 3.5 V Rechargeable Handle (wall plug-in)



#71630



#71670

#71000-C

#71000-A

3.5 V DESK CHARGER AND HANDLES

- #71140 (#71145* export) Universal Charger only
- #71670 Rechargeable 3.5 V Handle only
- #71410 Optional Wall Bracket for use with above desk chargers
- #72200 NiCad Rechargeable Battery (color-coded black)

3.5 V CONVERTIBLE HANDLE



- #71000-C (#71055-C* export) 3.5 V Convertible Power Handle (C-Cell battery capability)
- #72300 NiCad Rechargeable Battery (color-coded orange)

3.5 V Handle



- #71000-A (#71055* export) Rechargeable 3.5 V Power Handle
- #71050 Extra bottom section only
- #72300 NiCad Rechargeable Battery (color-coded orange)

*Available for export. Please specify country and voltage. Export ordering example: #71105 – Spain – 220 V

Diagnostic Sets

Welch Allyn offers a variety of customized equipment options to help fit your specific needs.

- Rechargeable Power Handle(s) of Choice
- Coaxial-Plus Ophthalmoscope and/or PanOptic™ Ophthalmoscope
- Halogen HPX™ Streak Retinoscope
- Carrying Case
- Fixation Card Set



#18335-SM Smart Set



#18335-C Classic Set



#18320-PC PanOptic Classic Set

#18335-SM Smart Set with Elite Retinoscope, Coaxial-Plus Ophthalmoscope, 3.5 V Lithium-Ion handle, and A/C charging module

#18335-C Classic Set with Elite Retinoscope, Coaxial-Plus Ophthalmoscope, and NiCad rechargeable convertible handle

#18320-PC PanOptic Classic Set with retinoscope, PanOptic Ophthalmoscope, transilluminator, and 3.5 V NiCad rechargeable convertible handle

	RETINOSCOPES		FIXATION CARDS	OPHTHALMOSCOPES				POWER HANDLES	CASES	
DIAGNOSTIC SET NUMBERS	#18200 SERIES STREAK	#18300 SPOT	#18250	#11735	#11730	#11720	#11820	#71900	#71000 or #71000-C**	HARD
#18335-SM	■		■	■				■		■
#18335	■		■	■					■	■
#18330-SM	■		■		■			■		■
#18330	■		■		■				■	■
#18430		■	■		■				■	■
#18320-PSM	■		■				■	■		■
#18320-PC	■		■				■		■	■
#18320-SM	■		■			■		■		■
#18320	■		■			■			■	■
#18420		■	■			■			■	■
#11790					■				■	■
#11770						■			■	■
#18340-SM	■		■					■		■
#18440		■	■						■	■

**Available with convertible handle, add "C" to set number when ordering (i.e. 18335 to 18335-C). All sets available in soft case. If soft case is desired, specify by changing last numerical digit of set number to a "1" (e.g., 18335 to 18331).



#05829-U Hard Case



#05620-U Soft Case



#05214-U Soft Case

767 Wall Transformer

The 767 Wall Transformer is designed to put everything you need for fast, accurate diagnoses at arm's length—and in one convenient modular package. The wall transformer holds two Welch Allyn handles, and features OptiSense™—a reliable optical sensor that automatically turns handles on when they're removed from the cradle and off when they are replaced. Best of all, you can add any Welch Allyn 3.5 V instrument to your system by simply attaching the optional third handle module.


- Handle-based rheostats allow for quick light-intensity adjustments
- Comfortable, contoured handles are easy to hold and maneuver
- Modular coiled cords are guaranteed for 10 years
- High-impact ABS housing is durable and easy to clean
- Improved cradle design featuring OptiSense

#76720 (#76725* export) 3.5 V Wall Transformer with clock (heads sold separately)

#76710 (#76715* export) 3.5 V Wall Transformer without clock (heads sold separately)

#76730 Third Handle Module

*Available for export. Please specify country and voltage. Export ordering example: #76725 – Italy – 220 V



Binocular Indirect Ophthalmoscopes

The Welch Allyn Binocular Indirect Ophthalmoscope (BIO) is one of the most advanced BIOs available. It offers state-of-the-art, convenient features along with a lightweight, comfortable fit and increased stability.

- #12500 BIO with Cobalt Blue Filter (optical portion and headband)
- #12500-D BIO with Cobalt Blue Filter and Diffuser Filter (optical portion and headband)
- #12500-Y BIO with Yellow Filter (optical portion and headband)
- #12500-DY BIO with Diffuser Filter and Yellow Filter (optical portion and headband)
- #74350 (#74355* export) Wall/Desk Power Source mounts to any convenient surface
- #74360 (#74365* export) Portable Power Source with up to 60 minutes of on-time—can be used with recharging transformer for extended use
- #05120-U Large Carrying Case
- #05122-U Small Carrying Case
- #12550 Teaching Mirror
- #12060 Scleral Depressor
- #50120 Fundus Charts and Pencils
- #73300 4 ft Cord Extension
- #12551 Voltage Adapter, 4.5 V

SMALL PUPIL CAPABILITY:
Convergence control allows illumination and visual pathways to be synchronized for excellent viewing in patients with pupils as small as 2 mm.




#74350 Wall/Desk Power Source



#74360 Portable Power Source

10-YEAR WARRANTY ON OPTICAL COMPONENTS AND HEADBAND

PADDED HEADBAND:

Well-padded leather contact points and "swooped" design provide maximum comfort and viewing stability. Easy to adjust for a secure fit.

ERGONOMICALLY DESIGNED CONTROLS:

Convenient, durable design provides simple operation and protection against accidental drops



#05120-U Carrying Case



SOFTILT PARK/REST FEATURE:

Allows optics to be moved to a stationary position out of your line of vision and eliminates risk of damage to eyeglasses or optics when returning them to use

HALOGEN LIGHT SOURCE:

One of the brightest, longest-lasting lamps available in a BIO (average lamp life 100 hours)

FULL-TIME SAFETY FILTER:

Protects your patients against UV and IR hazards

MIRROR TILT CONTROL:

Optimally synchronizes illumination and viewing paths

SEALED OPTICS:

Vital optics and mirror assemblies are sealed against dirt and dust for enhanced viewing and durability

IPD CONTROL:

Adjusts to handle interpupillary distances of 48 mm to 75 mm

OPTIONAL DIFFUSER FILTER:

Broadens the illumination beam for enhanced viewing of the peripheral retina

PRECISE ALUMINUM CHASSIS:

Protects optics and keeps components in precise alignment

SureSight™ Autorefractor

Portable, fast and affordable device that makes testing even your toughest patients quick and painless - for both you and your patients.

- Five-second, automatic testing is both fast and efficient
- Lights and sounds engage the patient's attention
- Measures and indicates abnormal readings automatically—no provider interpretation is required
- Minimal cooperation is required
- The hand-held, 2-pound unit and printer can be carried easily from room to room
- Easy to learn and use
- Refracts babies, children, and adults
- Test patients while they're wearing glasses or contact lenses

#14000 SureSight Vision Screener

#S1-14000-1 One-Year Comprehensive Partnership Program

#S1-14000-2 Two-Year Comprehensive Partnership Program



#14000



#43300

#41100

3.5 V Halogen Fiber Optic Transilluminators

General fiber-optic examination light provides cool illumination that eliminates patient discomfort. Useful in checking pupil reflex.

- Fiber optic light transmission for cool distal light
- Fully sterilizable with steam, water, gas or solutions
- Optional cobalt blue filter

#41100 Finnoff Ocular Transilluminator

#41101 Finnoff Ocular Transilluminator with Cobalt Blue Filter

#41102 Cobalt Blue Filter Only

#43300 Curved All-purpose Transilluminator

#03100-U Halogen HPX™ Replacement Lamp

#03100-LED LED Replacement Lamp



#76600

Halogen Professional PenLite

A convenient pocket light designed to fit your hand perfectly. Features bright, white halogen illumination for accurate tissue color rendition. Brass with black matte finish.

#76600 Halogen Professional PenLite

#07600-U Halogen Replacement Lamp



#34XFHT-B

Connex® ProBP™ 3400 Digital Blood Pressure Device

Accurate, reliable, connected blood pressure readings in the palm of your hand.

- Small, handheld design makes device easy to carry and capture readings in a triage setting
- Hospital-grade SureBP™ technology enables fast, accurate 15-second blood pressure readings, improving patient comfort
- Obtains accurate readings and minimizes error codes, even in the presence of many cases of weak pulse
- USB-wired technology enables connectivity to select EHR/EMR/HIS systems, improving documentation accuracy and efficiency
- Handheld, wall mount and mobile stand configurations
- Compatible with FlexiPort® blood pressure cuff sizes ranging from small child (size 8) to thigh (size 13)
- Lithium ion battery technology delivers 100 readings on a single charge with rapid recharge time
- Upgradeable so that your initial investment is protected as your needs change

Handheld

- #34XFHT-B Connex ProBP 3400 with SureBP
- #34XXHT-B Connex ProBP 3400 with Standard BP

Wall Mount

- #34XFWT-B Connex ProBP 3400 with SureBP
- #34XXWT-B Connex ProBP 3400 with Standard BP

Mobile Stand

- #34XFST-B Connex ProBP 3400 with SureBP
- #34XXST-B Connex ProBP 3400 with Standard BP
- #S1-3400 One-Year Comprehensive Partnership Program
- #S1-3400-2 Two-Year Comprehensive Partnership Program



#34XFTWT-B



#34XFST-B

Spot Vital Signs® Device

The Welch Allyn Spot Vital Signs takes vital signs at an affordable price. The fully automated, multiparameter Spot Vital Signs has become the device of choice for frontline caregivers.

- Fast and accurate—completes NIBP, pulse rate, temperature, and SpO₂ in approximately 30 seconds
- Multiparameter and modular; select the vital signs you need
- Connectivity to integrate with electronic health records (EHR)
- Set it on a countertop, mount it on a wall or wheel it on a mobile stand
- Big, bright, easy-to-read LCD screen
- One of the lowest-priced multiparameter devices on the market
- Fully upgradeable to meet your changing needs

#4200B-E1 Spot with NIBP only

#420TB-E1 Spot with NIBP and Temperature

#42M0B-E1 Spot with NIBP and SpO₂ (Masimo®)

#42MTB-E1 Spot with NIBP, SpO₂ (Masimo®) and Temperature

#42N0B-E1 Spot with NIBP and SpO₂ (Nellcor®)

#42NTB-E1 Spot with NIBP, SpO₂ (Nellcor®) and Temperature

#4700-60 Mobile Stand with Basket

#4701-62 Wall Mount with Basket



#42MTB-E1

#4200-170 Connectivity Accessory Kit with PC IR Adaptor

#05031-101 Disposable SureTemp Plus Thermometer Probe Covers (1,000 each)

#S1-4200 One-Year Comprehensive Partnership Program

#S1-4200-2 Two-Year Comprehensive Partnership Program

Spot Vital Signs® LXi Device

The Spot Vital Signs LXi is a faster, smarter, more efficient way to take vital signs. No matter how you configure it, the Spot Vital Signs LXi lets you gather the data you want with a clear, logical interface.

- SureBP™ measures blood pressure in about 15 seconds as the cuff is inflating—resulting in time savings for you and improved comfort for your patients. Pulse rate and MAP included
- Choice of SureTemp® Plus or Braun ThermoScan® PRO 4000 ear thermometer
- Optional Masimo® or Nellcor® SpO₂
- 50-reading memory allows retrieval of previous data for review or to download later
- Comes connectivity-ready, allowing you to connect with EHR systems
- Bigger, brighter LCD display provides greater visibility from farther distances
- Print your results or send them electronically

#450T0-E1 Spot LXi with SureBP and SureTemp Plus Thermometer

#450E0-E1 Spot LXi with SureBP and Braun ThermoScan® PRO 4000 Thermometer

#45MT0-E1 Spot LXi with SureBP, SpO₂ (Masimo) and SureTemp Plus Thermometer

#45ME0-E1 Spot LXi with SureBP, SpO₂ (Masimo) and Braun ThermoScan® PRO 4000 Thermometer

#45NT0-E1 Spot LXi with SureBP, SpO₂ (Nellcor) and SureTemp Plus Thermometer

#45ME0-E1 Spot LXi with SureBP, SpO₂ (Nellcor) and Braun ThermoScan® PRO 4000 Thermometer

#4700-60 Mobile Stand with Basket

#4701-62 Wall Mount with Basket

#4500-500 External Printer

#4500-925 Cable for Wired Connectivity

#05035-800 Braun ThermoScan® Disposable Thermometer Probe Covers (800 each)

#S1-4500 One-Year Comprehensive Partnership Program

#S1-4500-2 Two-Year Comprehensive Partnership Program



#4700-60

#45MT0-E1

Connex® Vital Signs Monitor 6000 Series

The Connex VSM is an advanced monitor that features a bright touchscreen display that is visually intuitive to operate and provides quick access to all device settings and controls—perfect for surgery centers and other ambulatory care offices. Choose the specific workflow you need—for both continuous monitoring and spot checking of neonatal, pediatric and adult patients.

- 15-second SureBP™ NIBP
- Optional SureTemp® Plus or Braun ThermoScan® PRO 4000 thermometry
- Optional Masimo® or Nellcor® pulse oximetry (SpO₂)
- Blood pressure intervals can be set for predefined or customized cycle times
- Only one step away from the home screen for any function or parameter
- Alarm controls are displayed and accessible on the home screen
- Save vitals information with the patient name for future reference
- Operate from an AC outlet or a high-performance lithium-ion battery
- 100% recharged in only four hours
- Connects to many leading EHR systems to reduce documentation errors
- Newly designed accessory cable management stand



Connex VSM 6300—1 host USB port for accessories, 6-Cell Li-Ion Battery*

- #63NXXE-B CVSM 6300 SpO₂ (Nellcor), Temperature (Braun ThermoScan PRO 4000)
- #63MXXE-B CVSM 6300 SpO₂ (Masimo), Temperature (Braun ThermoScan PRO 4000)
- #63NXPE-B CVSM 6300 SpO₂ (Nellcor), Temperature (Braun ThermoScan PRO 4000), Printer
- #63MXPE-B CVSM 6300 SpO₂ (Masimo), Temperature (Braun ThermoScan PRO 4000), Printer
- #63NTXX-B CVSM 6300 SpO₂ (Nellcor), Temperature (SureTemp Plus)
- #63MTXX-B CVSM 6300 SpO₂ (Masimo), Temperature (SureTemp Plus)
- #63NTPX-B CVSM 6300 SpO₂ (Nellcor), Temperature (SureTemp Plus), Printer
- #63MTPX-B CVSM 6300 SpO₂ (Masimo), Temperature (SureTemp Plus), Printer

Connex VSM 6400—4 host USB port for accessories, 9-Cell Li-Ion Battery*

- #64NXXE-B CVSM 6400 SpO₂ (Nellcor), Temperature (Braun ThermoScan PRO 4000)
- #64MXXE-B CVSM 6400 SpO₂ (Masimo), Temperature (Braun ThermoScan PRO 4000)
- #64NXPE-B CVSM 6400 SpO₂ (Nellcor), Temperature (Braun ThermoScan PRO 4000), Printer
- #64MXPE-B CVSM 6400 SpO₂ (Masimo), Temperature (Braun ThermoScan PRO 4000), Printer
- #64NTXX-B CVSM 6400 SpO₂ (Nellcor), Temperature (SureTemp Plus)
- #64MTXX-B CVSM 6400 SpO₂ (Masimo), Temperature (SureTemp Plus)
- #64NTPX-B CVSM 6400 SpO₂ (Nellcor), Temperature (SureTemp Plus), Printer
- #64MTPX-B CVSM 6400 SpO₂ (Masimo), Temperature (SureTemp Plus), Printer

*Other model configurations are available.

Connex Accessories

- #4800-60 Accessory Cable Management Stand
- #4701-62 Wall Mount with Basket
- #4500-925 USB Cable for wired connectivity
- #S1-6000 One-Year Comprehensive Partnership Program
- #S1-6000-2 Two-Year Comprehensive Partnership Program



Propaq® CS Monitor

The Welch Allyn Propaq CS is a cost-effective continuous vital signs monitor that features a bright, full-color touchscreen display, easy-to-use menu functions and is ideal for procedural sedation and recovery.

- 3- or 5-lead ECG, Masimo SET® or Nellcor® OxiMax® pulse oximetry, NIBP with Welch Allyn FlexiPort® technology, heart/pulse rate, respiration, temperature, capnography (sidestream and/or mainstream), invasive blood pressure (one or two channels)
- Full patient alarms including apnea alarming
- Trending on all parameters
- Programmable defaults on displayed parameters and alarm limits
- Printer

Please note that some of the features listed above are optional and should be discussed with your Welch Allyn representative.

#9001-001288 ECG, Respiration, Nellcor SpO₂, NIBP, 2 Temperature Channels, Battery, Connectivity Port

#9001-004720 ECG, Respiration, Nellcor SpO₂, 2 Temperature Channels, Printer, battery, Connectivity Port

#9001-005578 ECG, Respiration, Nellcor SpO₂, Sidestream CO₂, NIBP, 2 Temperature Channels, Printer, Battery, Connectivity Port

#S1-24X One-Year Comprehensive Partnership Program

#S1-24X-2 Two-Year Comprehensive Partnership Program



#802LT0N-UE1

Propaq® LT Monitor

Weighing in at just about two pounds, and with shock resistance that can withstand a six-foot drop, the durable Propaq LT is ideal for reliable spot-check and continuous monitoring needs during pre-op, post-op and simple procedures.

- 3- or 5-lead ECG, Masimo SET® or Nellcor® OxiMax® pulse oximetry, NIBP with Welch Allyn FlexiPort™ technology, heart/pulse, respiration
- Continuous monitoring and spot check modes
- Storage of up to 20 full-disclosure data snapshots
- Print 8.5" x 11" reports from your PC via USB
- Ability to attach PDF to an electronic medical record from your PC
- 24-hour tabular and graphical trending
- Compact design, around the size of a paperback book
- Large color display option
- Numerous mounting options to match your facility workflow

#802LT0N-UE1 ECG, Respiration, Nellcor SpO₂, NIBP, Charging Cradle with USB, Battery

#802LT0S-UE1 ECG, Respiration, Masimo SpO₂, NIBP, Charging Cradle with USB, Battery

#S1-802LT One-Year Comprehensive Partnership Program

#S1-802LT-2 Two-Year Comprehensive Partnership Program



#008-0931-00

Tycos® DS58 Classic Aneroid

Unsurpassed reliability in a traditional, handheld design with gear-free DuraShock™ technology.

- Certified accuracy to ± 3 mmHg
- Meets twice the AAMI shock-resistant standard for enhanced durability (60")
- Thumbscrew bleed valve
- Laser-engraved dial face for accuracy
- Jeweled movement for long life
- Premium inflation system for comfort
- Latex-free for safety
- Lifetime calibration warranty



#DS58-11
(Adult cuff
not shown)

- #DS58-11 Gauge with Reusable One-Piece Adult Cuff and Zippered Case
- #DS58-11CB Gauge with Two-Piece Adult Cuff and Bladder and Zippered Case
- #DS58-MC Multi-Cuff Kit with Reusable One-Piece Cuff (child print, small adult, adult, large adult) and Zippered Case
- #DS58-MCCB Multi-Cuff Kit with Two-Piece Cuff and Bladder (child print, adult, large adult) and Zippered Case

DS66 Trigger Aneroid

Innovative trigger-release design for easy operation.

- Certified accuracy to ± 3 mmHg
- Meets AAMI, shock-resistant standard for enhanced durability (30")
- Laser-engraved dial face for accuracy
- Jeweled movement for long life
- Trigger deflation valve for smooth, one-finger operation
- Latex-free for safety
- 15-year calibration warranty



#5098-27

- #5098-27 Gauge with Reusable One-Piece Adult Cuff
- #5098-27 CB Gauge with Two-Piece Adult Cuff and Bladder
- #5098-30 Multi-Cuff Kit with Reusable One-Piece Cuff (child print, small adult, adult, large adult) and Zippered Case
- #5098-30 CB Multi-Cuff Kit with Two-Piece Cuff and Bladder (child print, adult, large adult) and Zippered Case



#DS45-11

DS45 Integrated Aneroid

The first gear-free, shock-resistant aneroid sphygmomanometer with integrated cuff.

- Certified accuracy to ± 3 mmHg
- Meets AAMI shock-resistant standard for enhanced durability
- Laser-engraved dial face for accuracy
- Gauge rotates 360° for easy viewing and snaps directly into cuff for quick cuff change
- Latex-free for safety
- 10-year calibration warranty

- #DS45-11 Aneroid with Reusable One-Piece Adult Cuff
- #DS45-11C Aneroid with Reusable One-Piece Adult Cuff and Zippered Case
- #DS45-11CB Aneroid with Reusable Two-Piece Adult Cuff and Bladder



Advanced DuraShock™ Technology

Welch Allyn DuraShock Technology encompasses a gear-free design that is lighter and has fewer moving parts than traditional gauges. This unique design allows gauges with DuraShock Technology to be shock-resistant, providing additional durability and requiring calibration less often than many competing gauges. In fact, Welch Allyn gauges can fall 30 inches or 60 inches (depending on the gauge) onto a hard surface and still remain accurate. These gauges are perfect for the rugged, fast-paced office practice, and will provide more accurate readings with lower lifetime cost than traditional technology.

WELCH ALLYN OFFERS MERCURY DISPOSAL. Mercury spills are a costly and hazardous distraction. In fact, some parts of the country are regulating and even banning medical products containing mercury, such as thermometers and sphygmomanometers. Welch Allyn is helping—with our own mercury-free sphygmomanometers, and with a program that will safely dispose of your existing mercury sphygmomanometers. For more information, contact your local Welch Allyn representative or visit the EPA website at www.purdue.edu/dp/envirosoft/mercury/src/title.htm



Professional Adult Stethoscopes

Lightweight, comfortable, and affordable.

- General-use stethoscope
- Stainless steel double-head (bell and flat diaphragm) chestpiece with nonchill rim
- Rotatable binaurals and interchangeable comfort sealing ear tips
- 5-year warranty
- Late-free

- #5079-135 Single Lumen Professional Stethoscope (Adult, Double Head, Black)
- #5079-137 Single Lumen Professional Stethoscope (Adult, Double Head, Burgundy)
- #5079-139 Single Lumen Professional Stethoscope (Adult, Double Head, Navy)
- #5079-285 Single Lumen Professional Stethoscope (Adult, Double Head, Forest Green)
- #5079-145 Single Lumen Professional Stethoscope (Pediatric, Double Head, Black)



Lightweight Double-Head Stethoscope

This lightweight, high-quality stethoscope is perfect for blood pressure measurements and general auscultation. Featuring flexible single-channel tubing, it ensures extreme comfort even when used over extended periods

- Adjustable chrome-plated binaurals
- Lightweight and comfortable
- 31" overall length
- 2-year warranty

- #5079-73 Single Lumen Professional Stethoscope (Adult, Double Head, Velvet Black)
- #5079-74 Single Lumen Professional Stethoscope (Adult, Double Head, Poppy Red)
- #5079-75 Single Lumen Professional Stethoscope (Adult, Double Head, Dawn Gray)
- #5079-76 Single Lumen Professional Stethoscope (Adult, Double Head, Misty Blue)

OSZ-5 Automatic Self-Measurement Blood Pressure System

Take accurate blood pressure readings at home with this easy-to-use, automatic device.

- Accurate, easy-to-use automatic or manual blood pressure system
- Simple one-button operation is easy to use
- "Fuzzy logic" smart inflation system is fast and comfortable for the user
- Large, easy-to-read display shows systolic, diastolic, and pulse measurements simultaneously
- Optional voltage adapter to preserve batteries

- #7052-34 OSZ-5 Automatic Self-Measurement Blood Pressure System with Adult Cuff
- #7052-37 Voltage Adapter



Green Series™ Medical Exam Lights

Welch Allyn Green Series Medical Exam Lights are among the first to feature white LEDs that provide a brighter, whiter light than ever before and consume less energy. So now you can see truer tissue color during an exam or procedure while also reducing your current cost of ownership.

- 5,500 °K white light-emitting diodes (LEDs) as opposed to halogen, for a superior whiter, brighter light
- Truer tissue rendition and uniform spot—no dark or hot spots—for enhanced visualization of the exam area
- No bulb replacement for cost savings and less waste
- Consumes less energy for reduced cost of ownership
- Touchless auto on/off controls to reduce the risk of cross-contamination
- Ergonomic design that is easy to clean
- 5-year limited warranty



LED VS. HALOGEN

GS Exam Light IV

- Perfect for the OB/GYN office or for specialties that require very focused light in specific areas
- 3x the typical light output and features a highly flexible light pipe that is twice as long as conventional pipes
- Compact, unobtrusive head that ensures pinpoint accuracy and truer tissue color during an exam

- #48810 GS Exam Light IV with Table/Wall Mount
- #48950 Mobile Stand for GS Exam Light IV
- #48955 Table/Wall Mount
- #52640 Disposable Sheath (25/box, 5 boxes/case)



#44410

GS 300 General Exam Light

- Ideal for primary care facilities, outpatient facilities, dermatology and physician's office
- Perfect combination of intense light output and cool operation, all delivered in a compact design
- Affordable price and rugged performance

- #44400 GS 300 General Exam Light with Mobile Stand
- #44410 GS 300 General Exam Light with Table/Wall Mount
- #48955 Table/Wall Mount



GS 900 Procedure Light

- Perfect for minor procedures in the ambulatory care center, outpatient facility or physician's office
- Features 6 LEDs, outstanding durability, extraordinary maneuverability and unparalleled spot quality
- Ceiling mount option gives physicians the ability to maximize floor space and direct light wherever it is needed during a procedure
- Touchless intensity controls with minimal heat dissipation and high color temperature

- #44900 GS 900 Procedure Light with Mobile Stand
- #44900-C GS 900 Procedure Light with Ceiling Mount
- #44900-W GS 900 Procedure Light with Wall Mount
- #52630 Disposable Sheath (50/box, 5 boxes/case)


























GS 600 Minor Procedure Light

- Ideal for minor procedures and exams
- Features 3 LEDs and a wider head specifically for demanding viewing conditions that require exceptionally good light

- #44600 GS 600 Minor Procedure Light with Mobile Stand
- #44610 GS 600 Minor Procedure Light with Table/Wall Mount
- #48955 Table/Wall Mount
- #52630 Disposable Sheath (50/box, 5 boxes/case)




















Halogen Replacement Lamps

LAMP NUMBER	REPLACEMENT LAMP ACTUAL SIZE	CURRENT INSTRUMENTS	OBSOLETE INSTRUMENTS
00200-U		Transilluminators 41000 	
00900-U			Ophthalmoscopes Retinoscopes 10600 121 110 17600 17610 
01200-U		Binocular Indirect Ophthalmoscopes 125 Series 	
01300-U			PenLites 76000 77700 
01400-U			PenLite 77800 
01800-U			Retinoscopes 17200 
02300-U			Ophthalmoscopes 12900 
03000-U		Ophthalmoscopes StrabismoScope™ 11710 12400 	Ophthalmoscopes Retinoscopes 11600 11605 11610 18000 
03100-U		Illuminators Transilluminators 41100 43300 	
03400-U		Transilluminators 41110 	
03700-U			Retinoscopes 18100 

Halogen Replacement Lamps *continued*

LAMP NUMBER	REPLACEMENT LAMP ACTUAL SIZE	CURRENT INSTRUMENTS	OBSOLETE INSTRUMENTS
03800-U		PanOptic Ophthalmoscopes  11810 11820	
03900-U		Ophthalmoscopes  13010	Ophthalmoscopes  12810 11110
04400-U		Ophthalmoscopes  11470/11475	Ophthalmoscopes  11400 11411
04500-U			Retinoscopes  18010
04600-U			Retinoscopes  17710
04900-U		Ophthalmoscopes  11720 11730	Ophthalmoscopes  11620 11630
07200-U		Retinoscopes  17200	
07600-U		PenLite  76600	
08200-U		Streak Retinoscopes  18200 18245	Streak Retinoscopes  18235
08300-U		Spot Retinoscopes  18300	
08400-U		Retinoscopes  93300	

Replacement Batteries

BATTERY NUMBER	BATTERY	CURRENT INSTRUMENTS	OBSOLETE INSTRUMENTS
71960		Lithium-Ion Handle 71900 	
72300	3.5 V 	71000-A 71020-A 71000-C 71020-C 	71055-C 
72200	3.5 V 	71670 	71000 71020 
72000	2.5 V 	70700 71500 	60700 
72610		CompacSet™ 93400 	
72600		PocketScope™ 72800 	
72240		Binocular Indirect Ophthalmoscope 12000 74360 	

WELCH ALLYN LAMP POLICY

Welch Allyn instruments are completely engineered from the source of illumination to the instrument's usage. Lamp manufacture is a very critical part of this system; therefore, it is recommended that you use only genuine Welch Allyn lamps.

WELCH ALLYN BATTERY POLICY

Welch Allyn NiCad batteries will be replaced at a price that is reduced pro rata based on premature failure should they fail before the date printed on the battery (when used in Welch Allyn handles only). This date represents two years after the manufacture of the battery.

Welch Allyn Inc. Locations & Phone Numbers

NORTH AMERICA

Welch Allyn Inc.
Corporate Headquarters
4341 State Street Road
Skaneateles Falls, NY 13153-0220
USA

1 800 535 6663
Fax: 315 685 4091

Welch Allyn Inc.
Jordan Road Facility
4619 Jordan Road
Skaneateles Falls, NY 13153
USA
315 291 4301

Beaverton Office
Welch Allyn Protocol, Inc.
8500 S.W. Creekside Place
Beaverton, OR 97008
USA
503 530 7500
Fax: 503 526 4200

Canada Office
Welch Allyn Canada Ltd.
160 Matheson Boulevard East #2
Mississauga, Ontario L4Z 1V4
Canada
1 800 561 8797
Fax: 1 905 890 0008

LATIN AMERICA

Miami Office
10205 NW 19th Street, Suite 100
Doral, FL 33172
USA
305 669 9003
Fax: 305 669 8951

Argentina Office
Mariela Weigandt
Buenos Aires
Argentina
(011) 15 6206 9782
+54 9 11 6206 9782
+1 786 235 8178

Brazil Office
Alameda Tocantins 882
Barueri, Alphaville - SP
CEP - 06455-020
Brazil
4133-0800
Fax: 4133-0801

Chile Office
Hernando de Aguirre #162, Oficina 504
Providencia, Santiago
Chile
(56-2) 335-9134
Fax: (56-2) 335-6762

Colombia Office
Carrera 15A #122-60
Bogota, Colombia
(57-1) 481-7739
Fax: (57-1) 629-9234

Mexico City Office
Aristóteles 87, Polanco
Miguel Hidalgo, D.F. 11560
Mexico
52 555 283-2470
Fax: 52 555 283-2488

Mexico Office
Welch Allyn de Mexico S. de R.L. de C.V.
Calle Emilio Flores #2471-A
Colonia Canon de Padre
Tijuana, Baja California C.P. 22203
Mexico
1 800 535 6663

EUROPE, MIDDLE EAST & AFRICA

France Office
Welch Allyn France, SARL
Immeuble Le Galilée - 10 rue Fontaine
Rouge
Val d'Europe
France
77700 CHESSY
01 60 09 33 66
Fax: 01 60 09 67 97

Germany Office
Welch Allyn Service GmbH
Zollerstrasse
72417 Jungingen
Germany
0 7477 92710
Fax: 0 7477 927190

Iberia Office
Welch Allyn International Ventures Sales
Office
General Varela 35
planta 2ª, oficina 10-11
28020 Madrid Spain
+34 91 299 03 76
Fax: +34 91 425 26 79

Ireland Office
Welch Allyn Ltd.
IDA Business Park - Johnstown
Navan, County Meath
Ireland
+353 46 906 7790
Fax: 0 3 46 906 7754

Italy Office
Welch Allyn Italia
S.r.l.
+39 0269682425
Fax: +39 0269633697

Middle East & S. Europe Office
Zeppou 31 Street
16675, Glyfada
Athens, Greece
Fax: +30 210 9619892

Netherlands Office
Welch Allyn BENELUX
Staalweg 50
2612 KK Delft
The Netherlands
31 15 750 5000
Fax: 31 15 750 5050

Poland Office
Welch Allyn Poland & Baltic States
ul. Wilanowska 46
60-688 Poznan
Poland
+48 61 828 14 32

South Africa Office
Welch Allyn South Africa (Pty) Ltd.
394 Surrey Ave.
Ferndale, Randburg 2194
Johannesburg, South Africa
27 11 777-7555
Fax: 27 11 777-7556

Sweden Office
Welch Allyn
Sverige Aktiebolag Verkstadsgatan 12
363 39 Höganäs
Sweden
0114642330220
Fax: 0444642333210

U.K. Office
Welch Allyn (UK) Ltd.
Cublington Road - Aston Abbotts
Buckinghamshire HP22 4ND
England
0 207 365 6780
Fax: 0 207 365 9694

ASIA & PACIFIC RIM

Australia Office
Welch Allyn Australia (Pty) Ltd.
P.O. Box 132, The Metro Centre, Unit 5,
38-46 South Street
Rydalmere, NSW 2116
Australia
0 2 9638 3000
Fax: 0 2 9638 3500

China Office
Room 1011, Central Plaza
No. 227 Huang-Pi Bei Rd.
Huang Pu District
Shanghai, China
81 21 6327 9631
Fax: 81 21 6327 9632

India Office
No 41, 5th Cross Street, VGP Layout Part 2
Palavakkam, Chennai 600041, India
+ 91 44 24512917
Fax: +91 44 24512940

Japan Office
Welch Allyn Japan K.K.
6F Kinseisha Bldg., 3-15 Kanda-nishikicho
Chiyoda-ku, Tokyo 101-0054
Japan
81 3 3219 0071
Fax: 81 3 3219 0074

Malaysia Office
Welch Allyn Malaysia Sdn. Bhd.
10A, Jalan SS3/33
Taman Universiti
47300 Petaling Jaya
Malaysia
+603 78753341
Fax: +603 78754210

Singapore Office
Welch Allyn Singapore PTE LTD
438B Alexandra Rd,
#01-03 Alexandra Techno Park
Singapore 119968
+65 6419 8100
Fax: +65 6291 5780



Advancing Frontline Care™

welchallyn.com

800.535.6663

Local Rep: _____

Phone Number: _____

4341 State Street Road, P.O. Box 220, Skaneateles Falls, NY 13153-0220 USA
(p) 800.535.6663 (f) 315.685.3361 www.welchallyn.com

©2011 Welch Allyn Printed in USA SM1856 Rev D

Nellcor® is a registered trademark of Nellcor Puritan Bennett Inc. ThermoScan® is a registered trademark of Braun GmbH.
Masimo® is a registered trademark of Masimo Corporation.