

WelchAllyn®



Welch Allyn  
MacroView™ Otoscope




A bigger, better, sharper view

Enhancing visualization of the ear, helping you perform routine exams with **ease and confidence.**



Integrated Wall System



Diagnostic Desk Set




Diagnostic Set

### Better clarity and definition of landmarks

MacroView optics increase **magnification by 30%** versus a traditional otoscope, giving you a nearly full view of the tympanic membrane. Coupled with the widest focus range, MacroView offers a **bigger, better, sharper view.**

**2X** FIELD OF VIEW



Standard Otoscope

MacroView Otoscope


### Welch Allyn LED Lamps

LED technology provides the benefits of LED while operating at a color similar to halogen.




#### WITH THE FAMILIARITY OF HALOGEN

##### SureColor™ LED



Optimized color temperature, similar to what you are used to with halogen

##### Most Competitors' LED




Light too bright or too white can wash out critical details

#### LED IS ENERGY EFFICIENT

3X the on-time vs. halogen with Welch Allyn rechargeable handles

Lower operating temperature intended for greater patient/clinician comfort



HALOGEN

LED



**7-YEAR WARRANTY**

# Welch Allyn MacroView Otoscope



Proprietary **"Tip Grip"** helps ensure ear speculum is secured and ejects used ear speculum for easy disposal

Default focus position; ideal for most ear examinations

Adjustable focus wheel helps you accommodate for variable ear canal lengths and farsighted eyes, zooming in or out to fine tune view

**Optional throat/nasal illuminator** adds one more tool to your practice; easy to adjust on 3.5 V handle

3.5V Lithium Ion rechargeable handle (shown) provides 2.4X more power compared to our Nickel Cadmium rechargeable handle

## ACCESSORIES



**52432-U**  
Kleenspec® Disposable  
Ear Specula



**52700**  
Slotted Specula  
Instrumentation Tip



**23804**  
Otoscope Insufflation  
Bulb and Tube



**24330**  
Sofseal Specula Tip

## Welch Allyn Digital MacroView™ Otoscope

Supports the diagnosis with better communication  
to enhance patient compliance

- Store digital images in patient records to support documentation and diagnosis
- Offers the ability to view clear, live and still images of small details such as tubes and hearing aids
- An excellent training aid to support clinical education





## Configurations

<b>77791-MX</b>	Welch Allyn Green Series™ 777 Integrated Wall Diagnostic System including MacroView Otoscope, Coaxial Ophthalmoscope, KleenSpec Disposable Specula Dispenser, Wall Aneroid Sphygmomanometer and SureTemp® Plus 690 Electronic
<b>77792-M</b>	Welch Allyn Green Series 777 Integrated Wall Diagnostic System including MacroView Otoscope, Coaxial Ophthalmoscope, KleenSpec Disposable Specula Dispenser and Wall Aneroid Sphygmomanometer
<b>77792-MNOBP</b>	Welch Allyn Green Series 777 Integrated Wall Diagnostic System including MacroView Otoscope, Coaxial Ophthalmoscope and KleenSpec Disposable Specula Dispenser
<b>71641-MS</b>	Welch Allyn 3.5 V Halogen HPX™ Diagnostic Set including MacroView Otoscope, Coaxial Ophthalmoscope, two Lithium Ion Rechargeable Handles and Universal Desk Charger
<b>71641-M</b>	Welch Allyn 3.5 V Halogen HPX Diagnostic Set including MacroView Otoscope, Coaxial Ophthalmoscope, two Nickel Cadmium Rechargeable Handles and Universal Desk Charger
<b>23920</b>	Welch Allyn 3.5 V Halogen HPX Digital MacroView Fiber-Optic Otoscope; Mini-USB Cable (10.0 ft/3.0 m) and set of Reusable Ear Specula
<b>Accessories</b>	
<b>23804</b>	Insufflation Bulb
<b>52700</b>	4.00 mm KleenSpec Slotted Disposable Specula Instrumentation Tip
<b>52432-U</b>	2.75 mm Universal KleenSpec Pediatric Disposable Ear Specula
<b>52434-U</b>	4.25 mm Universal KleenSpec Disposable Ear Specula
<b>24320</b>	SofSeal Specula Tip for Disposable Ear Specula, small
<b>24330</b>	SofSeal Specula Tip for Disposable Ear Specula, medium

For more information on the MacroView Otoscope, please contact your Welch Allyn sales representative or visit [www.welchallyn.com](http://www.welchallyn.com).

Welch Allyn, Inc.  
4341 State Street Road  
Skaneateles Falls, NY 13153 USA  
(p) 800.535.6663 (f) 315.685.3361



[www.welchallyn.com](http://www.welchallyn.com)

**WelchAllyn®**