

Configurations Available

CONNEX VSM	6400 (6700)	6500 (6800)
Parameters		
SureBP NIBP	●	●
SureTemp Plus Thermometry	●	●
Braun PRO 4000 Ear Thermometry	●	●
Masimo/Nellcor SpO ₂	●	●
EarlySense Contact-Free Motion	+	+
Covidien Microstream EtCO ₂	+	+
Masimo RRa	+	+
Masimo SpHb	+	+
Communications		
5 USB Ports	●	●
Ethernet	●	●
802.11 a/b/g Wireless Radio	+	●
Accessories		
Barcode Scanner	●	●
Weight Scale	●	●
EMR Connectivity	●	●
Printer	●	●

● Included ● Optional + Upgradeable

ACCESSORIES

- 4800-60 Accessory Cable Management mobile stand w/ basket
- 4701-62 Wall mount with basket
- 008-0834-01 GCX wall mount with channel
- 04000-900 Braun PRO 4000 Dock with Thermometer
- 6000-NC Nurse call cable
- 6000-916HS Medical 2D barcode scanner with coiled USB
- 6000-920 Wireless radio kit—internal 802.11 a/b/g
- 4500-925 USB cable for wired connectivity
- 660-0138-00 Ethernet patch cable 5 ft.
- 6000-BT0 One-year service contract, NIBP and temperature
- 6000-BTS One-year service contract, NIBP, temperature, SpO₂
- 4900-60 Connex Accessory Cable Management Stand, Extended housing
- 4900-62 Wall Mount with Basket, Extended housing



For more information about the Connex Vital Signs Monitor, please contact your local Welch Allyn representative, call 1.800.535.6663 to speak with a customer service representative, or visit welchallyn.com/connex.

Welch Allyn Corporate Headquarters
4341 State Street Road, P.O. Box 220, Skaneateles Falls, NY 13153-0220 USA
(p) 800.535.6663 (f) 315.685.3361

WWW.WELCHALLYN.COM

© 2014 WELCH ALLYN SM4104 REV A

WelchAllyn

Advancing Frontline Care™



Welch Allyn Connex® Vital Signs Monitor 6000 Series

Accurate, automatic spot-check vitals and continuous monitoring sent to your EMR from the bedside for improved patient safety and clinical decisions.

The next generation in vital signs monitoring with easy-to-use color touch screen, advanced vital signs technology to monitor vitals, respiration and motion. Plus EMR integration to give clinicians immediate access to the information they need—where and when they need it.



Connect to the future of vital signs today

Designed specifically for busy floors like yours to increase efficiency

- > Captures both spot-check vitals and continuous monitoring
- > EMR integration for accurate electronic vitals documentation
- > Vivid color touch screen simplifies training and operation
- > Comprehensive point-of-care documentation including customizable data fields
- > Advanced Patient ID including two-way name confirmation
- > Patient data protection for HIPAA compliance
 - Protect patient identifiable information with customizable options, like abbreviations
 - Allow automated deletion of patient data after sending to EMR
 - Apply rules to check record completeness before saving
- > Complete error, easy-to-understand messages for easy troubleshooting

Industry-Leading Parameters for accurate patient information

- > SureBP® Blood Pressure
 - Adult, pediatric, neonatal blood pressure modes
 - Customizable blood pressure monitoring programs help nurses follow hospital protocol
 - 15-second readings for greater patient comfort
- > SureTemp® Plus and/or Braun ThermoScan® PRO 4000 ear thermometry (optional—choose one or both)
- > Masimo® or Nellcor® SpO₂ algorithms (optional)—both sensors and signal processing
- > Masimo noninvasive continuous monitoring of total hemoglobin (optional)
- > Manual measurements entry (weight, height, respiration rate, pain) and documentation of modifiers, like body position and O₂ therapy
- > Customizable documentation fields to fully document at the point-of-care


Upgradeable design to help protect your investment

Designed for the med/surg environment, with easy-to-use, configurable touchscreen display. And future-proof your purchase with these upgrade options designed to help keep patients safe and improve outcomes:

- > Braun ThermoScan PRO 4000 Ear Thermometer
- > Health o meter® Digital Weight Scale connectivity
- > Covidien® Microstream® Capnography (etCO₂)
- > Masimo® Acoustic Respiration (RRa)
- > Proven EMR connectivity
- > EarlySense® contact-free sensing
- > Masimo® SpHb noninvasive hemoglobin
- > Best-in-class wireless technology

Advanced accessory management to simplify your life

- > Optimal organization of accessories
- > Helps prevent tangling and sensor damage
- > Cleaning supply storage encourages proper hygiene
- > LED lighting helps in dark rooms without disturbing patients



Protect your investment with service programs and tools

- > Flexible Partners in CareSM programs cover all your service needs
- > Fully customizable, including programmable location identifier and quick USB configuration cloning
- > PartnerConnect™—ready for software-driven repairs and upgrades

