

# Connex® Spot Monitor Mounting and Mobility Options

Vital signs data within your reach—across care settings and workflows

Mobile stands support clinical workflows requiring a shared device across multiple rooms.



	Standard Basket Stand Simple portability with a small footprint	Work Surface and Accessory Management Stand Organized storage and cable management, plus a compact mobile work surface footprint	Accessory Power Management Stand Mobile work surface and integrated battery to power through a busy day
<b>Part number*</b>	7000-MS3	7000-MWS	7000-APM
<b>Height</b>	124.46 cm (49") total 96.52 cm (38") to top of basket handle	119.38 cm (47") total 99.7 cm (39.25") to work surface	119.38 cm (47") total 104.14 cm (41") to work surface
<b>Base Diameter</b>	58.42 cm (23")	58.42 cm (23")	58.42 cm (23")
<b>Work Surface Area</b>	N/A	29.51 cm (11.62") wide 28.56 cm (11.25") deep	29.21 cm (11.5") wide 25.4 cm (10") deep
<b>Accessory Basket</b>	31.75 cm (12.5") wide 20.96 cm (8.25") deep 18.42 cm (7.25") high	29.85 cm (11.75") wide 20.96 cm (8.25") deep 19.05 cm (7.5") high 4 accessory bays	29.85 cm (11.75") wide 20.96 cm (8.25") deep 19.05 cm (7.5") high 4 accessory bays
<b>Weight</b>	13.8 kg / 30.5 lb (fully loaded with device and accessories)	15.4 kg / 34 lb (fully loaded with device and accessories)	18.4 kg / 40 lb (fully loaded with device and accessories)
<b>Power Accessory Functionality</b>	N/A	N/A	9 cell Lithium Ion Battery (Field replaceable, integrated power supply, integrated work light, battery life indicator display, two USB ports)
<b>List price</b>	\$537.70 CAD	\$609.45 CAD	\$987.65 CAD

\*Connex Spot Monitor sold separately.

Wall mount options are ideal for care settings needing a dedicated device in every exam room.



#### Wall Mount with Basket

Easy-to-access storage and cable management

#### Integrated Wall System

Complete set of physical exam tools in one integrated package

<b>Part number(s)*</b>	7000-GCX	77794-MCSM**, 77794-2MCSM***
<b>Dimensions</b>	Wall Channel: 48.26 cm (19") Pivot Arm: 20.32 cm (8") Wall Mount Utility Basket: 32.39 cm (12.75") wide x 19.05 cm (7.5") deep x 14.99 cm (5.9") high	Height: 30.48 cm (12") Length: 86.36 cm (34") Depth*: 5.08 cm (2") *Depth dimension is for wall board only
<b>List Price</b>	\$432.60 CAD	77794-MCSM: \$1,746.05 CAD 77794-2MCSM: \$2,270.00 CAD

The Connex Spot Monitor uses a standard VESA mount. Aftermarket mounting options are available.



Learn more at [www.welchallyn.ca](http://www.welchallyn.ca).

\*Connex Spot Monitor sold separately.

\*\* 77794-MCSM model includes mounting bracket for Connex Spot Monitor, Coaxial Ophthalmoscope (#11720), MacroView™ Otoscope (#23810), KleenSpec® Disposable Specula dispenser, vital signs accessory basket, vital signs monitor power supply (#7000-PS) and wall mount panel.

\*\*\* 77794-2MCSM model includes mounting bracket for Connex Spot Monitor, PanOptic™ Ophthalmoscope with cobalt-blue filter and corneal viewing lens (#11820), MacroView Otoscope (#23810), KleenSpec Disposable Specula dispenser, vital signs accessory basket, vital signs monitor power supply (#7000-PS) and wall mount panel.