



Welch Allyn  
Connex® Spot Monitor

# Welch Allyn Connex<sup>®</sup> Spot Monitor

Accurate, automatic vitals sent to your EMR from the point of care for improved patient safety and clinical decisions.

The monitor to turn to for routine vital signs acquisition, interval monitoring, and custom scoring. Plus EMR integration to give clinicians immediate access to the information they need—where and when they need it.



## Improve patient outcomes and satisfaction

- > Reduce transcription errors and omissions by automatically transferring patient vital signs and other data directly into your EMR
- > Customizable fields to fully document at the point of care, eliminating the step of also entering some data in computers
- > Compatible with a wide range of cleaning products to reduce the risk of cross-contamination
- > Industry-leading parameter technologies to ensure accuracy and performance across a diverse range of patient populations

## Accessory Power Management Stand with WhisperDrive™ technology

- > Exclusive power management features to ensure performance even if the prior shift forgets to recharge the device
- > WhisperDrive technology provides smooth, quiet maneuverability to help minimize patient disturbance
- > LED lighting to help clinicians work in dark rooms without disturbing patients
- > Small work surface to improve ergonomics and ensure there's always a convenient place to document notes
- > Organization of accessories to help prevent tangling and sensor damage
- > Cleaning supply storage to encourage proper device cleaning

## Upgradeable design to help protect your investment

Designed for the med/surg environment with an easy-to-use, configurable touchscreen display. And future-proof your purchase with these upgrade options designed to help keep patients safe and improve outcomes:

- > Masimo<sup>®</sup> SpO<sub>2</sub>
- > Nellcor<sup>®</sup> OxiMax™ SpO<sub>2</sub>
- > Nonin<sup>®</sup> SpO<sub>2</sub>
- > Braun PRO 6000 thermometry
- > Barcode scanning
- > Compatible with a wide range of wireless networks and security protocols
- > Access to remote software upgrades and USB accessories




### Improve staff satisfaction and productivity

- > Designed by clinicians to increase efficiency and staff satisfaction
- > Intuitive touchscreen user interface to simplify training and operation
- > On-screen guides to prompt inexperienced users in navigating clinical workflow

### Secure patient information

- > Optional clinician sign-on screen restricts access to protected health information
- > Optional features, like abbreviations, help protect patient-identifiable information during clinical use
- > Automated deletion of patient data after sending to EMR



### Improve patient safety

- > Captures custom patient scores at the bedside to allow clinicians to document and track signs of patient deterioration
- > Allows the upload of up to three protocols the facility already uses—such as MEWS or other calculations
- > Displays facility's own recommendations based on patient scoring results to help guide clinicians through next steps and to make informed care decisions quickly



### Partners in Care<sup>SM</sup> Services for Connex Spot Monitor

Welch Allyn offers a full suite of services to help protect your organization's investment and ensure clinicians have what they need to make informed patient care decisions.

- > **Unique Enterprise Deployment Solution**—to ensure your monitors arrive configured to your specification right out of the box
- > **Biomedical Service and Repair Training**—to empower staff to maintain and repair devices in-house
- > **Support Programs**—cover repair parts and replacement accessories for efficient, cost-effective device support
- > **Service Dashboard**—enables proactive monitoring of connected devices, configuration creation, and remote deployment

## Configurations Available

CONNEX SPOT MONITOR	7100	7400	7500
<b>Parameters</b>			
SureBP® NIBP	●	●	●
SureTemp® Plus Thermometry	○	○ ↑	○ ↑
Braun PRO 6000 Ear Thermometry	○	○ ↑	○ ↑
Masimo or Nellcor OxiMax® SpO <sub>2</sub>		●	●
Nonin SpO <sub>2</sub>	○		○
<b>Communications</b>			
WiFi		↑	●
USB, Ethernet	●	●	●
<b>Accessories</b>			
Classic Stand	○	○	○
APM Stand	○	○	○
GCX Wall Channel	○	○	○
Barcode Scanner	○	○	○

● Included    ○ Optional    ↑ Upgradeable

### ACCESSORIES

<b>7000-APM</b>	Connex Spot Accessory Power Management Stand
<b>7000-MS3</b>	Connex Spot Classic Mobile Stand
<b>7000-GCX</b>	Connex Spot GCX VESA Wall Channel
<b>6000-NC</b>	VSM 6000 Nurse Call Cable
<b>7000-916HS</b>	HS1 2D Barcode Scanner
<b>4500-925</b>	Connectivity Accessory Kit USB Cable
<b>6000-50</b>	VSM 6000 USB Configuration Memory Stick
<b>7000-PS</b>	Power Supply (required if not ordering a stand or wall mount)
<b>BATT22</b>	CSM Lithium-Ion Battery
<b>BATT99</b>	APM Stand Lithium-Ion Battery
<b>4500-35</b>	Blood Pressure Hose with FlexiPort (10 ft)
<b>7000-33</b>	Neonatal Blood Pressure Hose (10 ft)
<b>LNCS-DCIP</b>	Masimo Pediatric Reusable Finger Sensor
<b>D-YS</b>	Nellcor Dura-Y® Sensor
<b>D-YSPD</b>	Nellcor PediCheck® Pediatric Finger Sensor (Requires D-YS)
<b>2360-010</b>	Nonin Reusable Pediatric Sensor (2 m)
<b>S1-CSM-5</b>	CSM, Comprehensive Partner Program, 5 years
<b>S2-CSM-5</b>	CSM, Biomed Partner Program, 5 years
<b>S4-CSM-5</b>	CSM, Warranty Extension, 5 years

*One- and two-year service plans also available.*

Welch Allyn, Inc.  
 4341 State Street Road  
 Skaneateles Falls, NY 13153 USA  
 (p) 800.535.6663 (f) 315.685.3361

[WWW.WELCHALLYN.COM](http://WWW.WELCHALLYN.COM)

**WelchAllyn®**  
 Advancing Frontline Care™