



# RIGID ENDOSCOPY

**WelchAllyn**<sup>®</sup>

Advancing Frontline Care<sup>™</sup>

# Quality and Value in Ri

Welch Allyn offers a wide range of rigid endoscopy products for anoscopy and sigmoidoscopy, complemented by versatile fiber optic illumination systems for optimum visualization.

The sigmoidoscope range includes both reusable and disposable instruments and a variety of sizes to accommodate differing requirements. Welch Allyn fiber optic anoscopes represent the finest in design and craftsmanship for optimal performance.



## Welch Allyn Disposable Sigmoidoscopes and Anoscopes

One of our most popular designs, the Welch Allyn Disposable Sigmoidoscope may be used in conjunction with either a 6 V illumination system or traditional Welch Allyn 3.5 V power handles.

- > Highly visible graduations to assist in positioning
- > Glare-free frosted coating ensures maximum light output at distal tip
- > Individually wrapped for safety and convenience
- > Fiber-optic illumination system provides a cool 360 degree ring of light
- > Hinged window with neoprene seal to avoid fogging during insufflation

<b>53130</b>	KleenSpec® Disposable Sigmoidoscope with Obturator (case of 100)
<b>36103</b>	Complete Illumination System for Disposables
<b>73210</b>	Light Handle with 0.9 m Cord
<b>73220</b>	Light Handle (73210) plus Transformer (73305)
<b>36019</b>	Fiber Optic Light Head
<b>73305</b>	Transformer with 1.5 m Cord
<b>73500</b>	3.5 V Halogen Handle Adapter
<b>71000-A</b>	3.5 V Direct Plug-in Rechargeable Handle
<b>30200</b>	Insufflation Bulb, complete
<b>30210</b>	Hydrophobic Filters (pack of 10) with Luer Lock

## Disposable Anoscope

<b>53110</b>	KleenSpec® Disposable Anoscope with Obturator (box of 25)
--------------	---

### Disposable Anoscope Accessories

<b>73210</b>	Light Handle with 0.9 m Cord
<b>73220</b>	Light Handle (73210) plus Transformer (73305)
<b>36019</b>	Fiber Optic Light Head
<b>73305</b>	Transformer with 1.5 m Cord



# gic Endoscopy



35303  
Fiber Optic Sigmoidoscope Set

## Standard and Pediatric/Stricture Sigmoidoscopes and Anoscopes

All Welch Allyn reusable sigmoidoscopes/anoscopes feature:

- > Hinged window with a neoprene seal to avoid fogging during insufflation
- > Graduated tube with removable tip for fiber cleaning (metal only)



### Standard Sigmoidoscopes

- 32810** 19 mm x 15 cm Standard Sigmoidoscope, Speculum only with Obturator
- 32820** 19 mm x 25 cm Standard Sigmoidoscope, Speculum only with Obturator
- 33830** 19 mm x 30 cm Standard Sigmoidoscope, Speculum only with Obturator
- 32830** 19 mm x 35 cm Standard Sigmoidoscope, Speculum only with Obturator
- 33220** 23 mm x 25 cm Standard Sigmoidoscope, Speculum only with Obturator



### Pediatric/Stricture Sigmoidoscopes

- 32010** 11 mm x 10 cm Fiber Optic/Stricture Speculum only with Obturator
- 32020** 11 mm x 25 cm Fiber Optic/Stricture Speculum only with Obturator
- 32410** 15 mm x 15 cm Fiber Optic/Pediatric Speculum only with Obturator
- 32420** 15 mm x 25 cm Pediatric/Stricture Speculum only with Obturator



### Rotating Anoscope

Anoscope Only 68 mm Speculum Length

- 38900** Rotating Speculum Anoscope with Obturator, 20 mm inside diameter, 22 mm outside diameter




### Fiber Optic Anoscope

Speculum Only with Obturator

- 37019** 19 mm x 7 cm Fiber Optic Anoscope, Speculum Only with Obturator
- 37023** 23 mm x 7 cm Fiber Optic Anoscope, Speculum Only with Obturator
- 37027** 27 mm x 7 cm Fiber Optic Anoscope, Speculum Only with Obturator

# Flexible Illumination Systems for Optimal Visualization



## Sigmoidoscope/Anoscope Lighting Systems

Powerful illumination systems.

- > Halogen light provides consistent, long-lasting illumination (73500)
- > Fiber optics provide a cool distal ring of light with no reflections or obstructions
- > Sigmoidoscopes/anoscopes can be used with a choice of illumination systems—a 6 V system with optional on/off switch and rheostat control, or with Welch Allyn's 3.5 V hand-held power sources (with adapters 73500 and 73550).
- > 73500 adapter fits most Welch Allyn 3.5 V power handles

## Sigmoidoscope/Anoscope Lighting Systems and Accessories

- 36103** Illumination System for Disposable Sigmoidoscopes (Includes: 73305, 73210, 36019 and 30200)
- 36019** Fiber Optic Light Head for Disposable Sigmoidoscopes and Anoscopes
- 38700** Fiber Optic Light Carrier for use with Standard Anoscopes and Light Adapter
- 73500** 3.5 V Halogen Handle Adapter

- 73305** Transformer with 1.5 m cord
- 73210** Light Handle with 0.9 m cord
- 03100-U** 3.5 V Halogen Lamp (for 73550)
- 30200** Insufflation Bulb, Complete
- 30130** 30 cm Suction Tube (Standard)
- 30140** 40 cm Suction Tube
- 07800-U** 6 V Halogen Lamp