



Welch Allyn  
HScripte® Holter Monitor Analysis System





## Freedom

### Empower your patients to go about their daily lives.

Ideal for traditional Holter monitoring, the H3+™ Digital Recorder provides 3-channel recordings using a 5-wire patient cable. Weighing in at one ounce (28 grams), the H3+ gives new meaning to the word compact, while providing some of the most impressive capabilities available today:

- Integrated waveform display
- Pacemaker detection
- 3-channel recording
- Up to 48 hours of recording using a single AAA battery





With the industry's broadest range of diagnostic cardiology solutions, we help people get better care, inside and outside the hospital. Backed by **clinical excellence, connected solutions and continuous innovation**, Welch Allyn Cardiology is proud to be powered by Mortara.



# Efficiency

Combine speed, flexibility and performance.



Data transfer and automatic analysis require **only about a minute**.

With intuitive menus and color-coded beat identifiers, the system quickly leads clinicians to the most significant arrhythmia events for review. Automatic strips enable rapid reporting of results, optimizing clinicians' time. Flexibility in scanning and report formatting streamlines clinical workflows and satisfies distinct user preferences. Scheduling services and Holter study tracking allow you to monitor study progression and manage workflow tasks.



**The best-in-class VERITAS™ ECG algorithms demonstrate impressive performance in beat detection and labeling.**

VERITAS provides accurate and meaningful results presented in trends, tables and histograms, to assist with determination of frequency and severity.



# Accessibility

## Information is powerful when it is accessible.

The HSCRIBE System eliminates the physical constraints of traditional Holter systems and brings forth the benefits of information sharing. No longer is the patient or recorder needed at your site. Instead, perform recorder hookup and download at your office, clinic or hospital and simply transfer the recordings over a secure internet or intranet connection.

The HSCRIBE database maintains centralized patient demographics and provides a repository for final reports, allowing clinicians to quickly search for patient information or a specific Holter record. Management of HSCRIBE system settings includes password-protected user permissions, application settings, audit trails for regulatory compliance and electronic signature, all of which are customizable based on facility preference.

Save time and reduce costs by sharing and expediting Holter study workflow and results between networked locations with bidirectional connectivity via industry-standard DICOM® or HL7 exchange.



Contact your Welch Allyn representative or visit [www.welchallyn.com](http://www.welchallyn.com) to learn more.



Welch Allyn, Inc.  
4341 State Street Road  
Skaneateles Falls, NY 13153 USA  
(p) 800.535.6663 (f) 315.685.3361