



# 313 Steel Frame Folding Screen Parts List



PARTS LIST				
Image Number	Note	Part Number	Description	QTY per Product
1	A	N/A	313 SCREEN FRAME	1
2		211650	SCREEN LEG W/CASTERS ASSEMBLY KIT (313-342)	2
2A	B	211610	DOUBLE LEG 313/342	2
2A	B	702220	CASTER INSERT, 1" RND X 16 - 18 GA	4
2A	B	400000	CASTER 2" BLACK/TWIN HOODED AP7	4
2B	B	566506	SCREW-MACH 1/4-20 X 1 1/2 RND PHIL M/S ZINC	4
2B	B	566204	NUT-1/4-20 HEX NYLON LOCKNUT ZINC	4
2C		701700	CRUTCH TIP 1" BLACK	4
NOTES				
Note	Description			
A	This part is not available for purchase			
B	Parts MUST be ordered in the above listed kit. These items can NOT be sold seperately.			

# 313 Steel Frame Folding Screen Parts List



PARTS LIST				
Image Number	Note	Part Number	Description	QTY per Product
1A		299800	313 VINYL PANEL STANDARD WHITE	3
1A		2998XX	313 SURE CHEK PANEL	3
1B		302005	313 SCREEN ROD, ALUMINUM	6
1C		704313	INSERT, SCREEN ROD PLASTIC, 313	12
1D		302015	SPRING, SCREEN ROD	6
NOTES				
Note	Description			
A	XX- For SURECHEK use color code for specific part number i.e. 93 - Mint (299893)			