

Designer Privacy Screen 3170



ISO 9001-2008 Certified

Designer quality privacy screen that will coordinate with every room.

Winco's Designer Privacy Screen has all the features you need, the reputation you trust, and a contemporary design to meet everyone's needs. Match the curtain color to any standard Winco vinyl color, and you have a custom designed privacy screen that coordinates with the other Winco products in your facility.



Use your own custom material to create the perfect atmosphere in your facility

KEY FEATURES

- Durable steel frame
- Twin casters
- Replaceable screen panels
- Hinged panels

Proudly Made in the USA



1.800.237.3377 • 5516 SW 1st Lane • Ocala, Florida 34474

www.wincomfg.com



PRODUCT SPECIFICATIONS

DESIGNER PRIVACY SCREEN 3170

- 1 Choose from our extensive vinyl collection including our Sure-Chek[®] antimicrobial material and Silvertex.
- 2 Hinged panel ensures years of problem-free use.
- 3 Wide base for added stability
- 4 Smooth rolling twin wheel casters with optional crush-tips included.
- 5 Smooth powder-coated steel frame for durability.



SPECIFICATIONS

3170

Base Size.	23" x 31" (58.4 x 78.7cm)
Overall Height.	70" (177.8cm)
Overall Width.	84" (213.4cm)
Single Panel Length.	55" (139.7cm)
Single Panel Width.	28" (71.1cm)
Weight.	45.5 lbs. (20.6kg)


COLORS

STANDARD COLORS



Printed colors shown may vary slightly from actual color of product.

OPTIONAL COLOR UPGRADES (ADDITIONAL CHARGES APPLY)



SILVERTEX VINYL (ADDITIONAL CHARGES APPLY)



SURE-CHEK[®] VINYL (STANDARD)



Sure-Chek[®] vinyls are antimicrobial, antistatic, odor resistant, stain resistant, and fire retardant.

01/31/11

Rev. B

009063