

Premier LifeCare Recliner 5400



®

ISO 9001-2008 Certified

The Winco 5400 is an upscale long-term recliner with fully upholstered sides for contemporary modern style, and practical function. Durable and comfortably contoured urethane armrests protect the upholstered finish and are easily cleaned. The 5400 is a perfect blend of value, function and style.



KEY FEATURES

- Urethane armrests
- Fully upholstered
- Multiple side-table Options
- 3-position lockbar
- "No pinch" safety zone

Proudly Made in the USA




1.800.237.3377 • 5516 SW 1st Lane • Ocala, Florida 34474

www.wincomfg.com



PRODUCT SPECIFICATIONS

PREMIER LIFE CARE RECLINER 5400



- 1 Push handle for easy maneuverability
- 2 Comfortable, contoured urethane armrests to keep upholstery free from damage.
- 3 Black, low profile 3" casters provide smooth and easy movement
- 4 3 Position Lock Bar feature allows for safe secure and comfortable positioning.



SPECIFICATIONS

5400


Seat Height	20.5" (52cm)
Seat Width	19.5" (49.5cm)
Seat Depth	20" (50.8cm)
Overall Height	46" (116.8cm)
Overall Width	27.75" (70.5cm)
Back Height Above Seat	30.5" (77.5cm)
Overall Length (Upright)	39.5" (100.3cm)
Overall Length (Reclined)	62" (157.5cm)
Weight	90 lbs. (40.8kg)
Weight Capacity	350 lbs. (158.7kg)

Add 1" to chair width for each side-table in down position
 *Add 10.5" to chair width for each raised side table

OPTIONS




STANDARD COLORS



Printed colors shown may vary slightly from actual color of product.

10-03-11

OPTIONAL COLOR UPGRADES (ADDITIONAL CHARGES APPLY)



SILVERTEX VINYL (ADDITIONAL CHARGES APPLY)



REV.E

009100

COLORS