



ISO 9001-2008 Certified

## Hassle-free design to optimize your patient's comfort

This Elite Care Cliner features dual swing-arms, allowing patients with limited mobility a modified entry. Also, LiquiCell® has been added into the seat for extra comfort. This ultra-thin, liquid-filled interface aids in improved blood flow, and helps prevent skin breakdown.



**HEAT  
&  
MASSAGE**  
(Optional)

## KEY FEATURES

- Dual side swing-arms for modified entry and ease of cleaning
- Dual fold down tables
- LiquiCell® for added comfort
- 5" nylon® total-lock casters
- Built-in headrest
- Headrest cover

Proudly Made in the USA



1.800.237.3377 • 5516 SW 1st Lane • Ocala, Florida 34474

[www.wincomfg.com](http://www.wincomfg.com)



# PRODUCT SPECIFICATIONS

## ELITE CARE CLINER STANDARD X-LARGE 6940/6950

- 1 Standard Trendelenberg release on back with easy-access handle
- 2 Wide elastic bands on back frame and seat provide enhanced comfort and support.
- 3 Non-porous, fold-down side tray contains small spills. STANDARD on right and left sides.
- 4 Dual Swing-arm armrests, w/ quick release, for easy entry and ease of cleaning
- 5 Nylon total-lock and directional-lock casters provide ease of maneuverability
- 6 Easy swing-arm release lever
- 7 **LIQUIGELL**® improves blood perfusion and aids in preventing skin breakdown.



**6950**  
XL Elite Care Cliner


SPECIFICATIONS

	<b>6940</b>	<b>6950</b>
Seat Height . . . . .	21" (53.3cm)	21" (53.3cm)
Seat Width . . . . .	21" (53.3cm)	25" (63.5cm)
Seat Depth . . . . .	20" (51cm)	20" (51cm)
Overall Height . . . . .	48.5" (123.2cm)	48.5" (123.2cm)
Overall Width . . . . .	33.5" (85.1cm) Add 9" for each side tray up	37.5" (95.3cm) Add 9" for each side tray up
Back Height Above Seat . . . . .	32.5" (82.6cm)	32.5" (82.6cm)
Overall Length (Upright) . . . . .	41" (104.1cm)	41" (104.1cm)
Overall Length (Reclined) . . . . .	76.5" (194.3cm)	76.5" (194.3cm)
Armrest Height . . . . .	28.5" (72.4cm)	28.5" (72.4cm)
Weight . . . . .	104 lbs. (47.2kg)	127 lbs. (57.6kg)
Weight Capacity . . . . .	350 lbs. (158.7kg)	450 lbs. (204kg)

OPTIONS


 <b>IB-133 IB</b>	 Heat <b>HT</b> Heat & Massage <b>HM</b>	 Armrest Covers <b>AC00</b>	 Side Cushion <b>SC00</b>	 IV Pole <b>IV</b> (Left rear only)	 Leg Extensions <b>SE00</b>
---	---	--	--	---	--

**STANDARD COLORS**



Printed colors shown may vary slightly from actual color of product.

**OPTIONAL COLOR UPGRADES**



**SILVERTEX VINYL**

(ADDITIONAL CHARGES APPLY)



COLORS