

Why the SUMMIT

What sets the Summit apart



Smallest Footprint on the Market

Telescoping Sling Bar

Lithium-Ion Battery Option

Fully Adjustable Shin Pad

O₂ Tank/IV Pole Holder

Removable Footplate

Adjustable Base

Encapsulated Casters



500 POUND
LOAD CAPACITY

Smallest Footprint on the Market: 23.8" width - fits in small places or departments without a lot of storage.

Telescoping Sling Bar: allows for easier sling application and operation with different patient body types.

Lithium-Ion Battery Option: averages 40 lifts per charge.

O₂ tank/IV Pole Holder: eliminates the need for a caregiver to follow with another device.

Fully Adjustable Shin Pad: adjusts up/down, in/out, or removable to accommodate various leg lengths.

Removable Footplate: the removable footplate allows users of 2 different mobility levels.

Adjustable Base: the Summit's base can spread for easy access to wheelchair, commode, and much more.

Encapsulated Casters: no hair or dirt buildup.

500lb Load Capacity: generous weight limit enables use by patients



800.255.3126 | wyeastmedical.com